

## PROCEDURE CHAIR

HT-100 DLX (3 actuator)  
/ HT-200 DLX (2 actuator)



## Applications

General examination, dermatology and cosmetology procedures, Laser treatments, Hair Transplantation.

## Specifications

- **Make & Model:** HT-100 DLX (3 actuator) / HT-200 DLX (2 actuator)
- **Material :** Epoxy Powder Coated Steel.
- Wired remote controlled Actuator based technology
- Body made of Powder coated mild steel
- Weight bearing capacity upto 150Kgs
- Height adjustable rotatable armrest
- Detachable handrests for sleeping position Tiltable height adjustable headrest
- PU foam for maximum patient comfort
- **Dimensions:** Length - 6Ft, Width-2 Ft height adjustable from 24 inches to 38inches
- This model can be fitted with either 2 or 3 actuators to control the movements
- One of the actuator controls up and down movement
- Second actuator controls synchronous movement of back rest to foot rest
- While the third actuators enable separate controls for backrest and footrest. The chair can be used for general examination, dermatological and cosmetologically procedures, Laser treatments, Hair Transplantation (**OPTIONAL**)
- **Product Dimension & Weight :** 72x24x24-38inch & 100Kg

NOTE : DUE TO OUR CONTINUOUS R & D PROGRAM SPECIFICATIONS AND DESIGN MAY CHANGE WITHOUT PRIOR NOTICE.



**C Cube Advanced Technologies**  
(An ISO-9001 & ISO-13485  
Certified Company)



+91-80500 50512



+91-80-43021028



marketing@ccubeadvtech.com  
ccubeadvtech@gmail.com



www.ccubeadvtech.com

Manufacturers, Suppliers & Dealers of Medical Equipments