

LASER HAIR GROWTH Laser Optimus

The product is a cutting edge technology with zero side effects and has delivered very good results. The product can be deployed at dermatology, skin care cosmetology centers. The product covers complete scalp promoting hair growth uniformly throughout scalp .

Features & Benefits

- Completely portable and easy to use.
- You don't have to hold the device (unlike a laser hair comb) Hands Free, completely fits in almost any hat.
- Get Thicker, Fuller, Healthier hair if you have thinning hair
- **Reverse hair loss** : effective for thinning hair and alopecia areata.
- Convenient and discreet.
- Use at home, while walking, driving anywhere!
- No more scheduled appointment for treatment.
- Complements other hair loss therapies.
- Hands free, convenient hair growth system provides full scalp coverage, is fully portable, easy to use and highly effective Low-Level Laser Therapy (LLLT) system that enables significant hair growth.
- Low Level Laser Therapy (LLLT) is a highly effective hair growth treatment process. LLLT increases cellular activity and promotes the proliferation of thick, full and healthy hair. The absorption of light delivered via an acute and measured wavelength (650nm) promotes hair growth.
- **GROWS HAIR** : Low-Level Laser Therapy (LLLT) technology utilizes highly effective laser diodes to effectively stimulate and energize cellular activity causing an uptake in the natural function of the hair follicle. This proven 4 to 6 month process reverses hair thinning and loss.
- **SAFE AND EFFECTIVE** : The Laser Optimus 28 / 60 is a safe and effective device, designed and produced specifically for the purpose of growing hair in people affected with Androgenetic Alopecia (thinning hair). Laser Optimus 28 is a non-invasive; Low Level Laser Therapy (LLLT) technology has no side effects, unlike pharmaceuticals, topicals and pills used to treat hair-loss.
- **AFFORDABLE** : Laser Optimus is not only highly effective. It's affordable as well. A onetime purchase will provide you with everything you need to grow and maintain strong, healthy hair for lifetime.
- **OPTIMUM SCALP COVERAGE** : The Laser Optimus 28 / 60 is designed and crafted to effectively cover the scalp and engineered to focus on areas prone to hair loss. The inner helmet canopy is designed to conform to contour and curvature of the scalp, ensuring consistent and equal distribution of laser for optimum results.



Safe to use
LLLT Technology



Advantages

- Hair Regrowth
- Hair Restore
- Prevent Hair fall
- Designed as per male & female hair loss pattern.

Specifications

- **Make & Model:** C CUBE, LASER OPTIMUS
- **Laser Diodes** : 28 Nos, 650nm 5mW, Light output through special laser multiplier.
- Cooling Fan provided for patient comfort
- **Main power supply** : 230 VAC /50Hz.
- **Battery** : Rechargeable Battery Operated with auto cut off feature after 21 minutes.
- **Body of the entire unit** : FRP, combination of cream and grey colours.
- Bag provided to keep and carry unit and to avoid the dust.
- User manual provided to operate the unit easily.
- **Warranty** : 6 months for the battery and one year for the unit.
- **Unit Dimension:** 29X23X20 cms
- **After packing** :34x27x25 cms
- **Unit Package weight** : 4 Kg.
- **Volumetric Weight** : 6 Kgs

Accessories

- Carry Bag
- Charging Adapter
- Pre Set Timer Controller
- User Manual

NOTE : DUE TO OUR CONTINUOUS R & D PROGRAM SPECIFICATIONS AND DESIGN MAY CHANGE WITHOUT PRIOR NOTICE.



C Cube Advanced Technologies
(An ISO-9001 & ISO-13485
Certified Company)

+91-80500 50512 +91-80-43021028

marketing@ccubeadvtech.com
ccubeadvtech@gmail.com

www.ccubeadvtech.com

Manufacturers, Suppliers & Dealers of Medical Equipments