

CENTRIFUGE MACHINE

A centrifuge is a system that separates the component parts of a liquid or fluid (and even gases) using centrifugal force. This is done by rotating the fluid in a bottle at high speed, thereby separating fluids of varying densities or liquids from solids.

These are low priced units meeting the majority of the small capacity centrifuging requirements of medical practitioners, pathological laboratories and hospitals for routine centrifuging tests. The units are fitted with resiliently mounted motor for vibration free performance having built in 5 step speed regulator. Available with Swing out or Angle rotor.

CF-100



- **Make & Model :** C CUBE, CF-100
- **Speed:** 4000 RPM
- **Speed accuracy:** +/-20
- **Time range:** 1~30 min.
- **Rotor capacity:** 6x15ml
- **Total volume:** 300ml
- **Type of head:** Swing out
- **No. of tubes:** 6
- **Tube size:** 15ml, 50ml
- **Max. RCF:** 1880 g
- **Power Supply :** 100 to 240.VAC, 50/60Hz, 70W, 3A, Single phase.
- **Dimension (LxHxW):** 280x266x364mm.
- **Weight:** 6 Kg.

CF-150



- **Make & Model :** C CUBE, CF-150
- Best grade material
- High performance
- Long functional life
- Reduced maintenance costs
- Constant performance
- Easy to operate
- Rectangular body with solid state Stepless speed regulator.
- Digital Speed (RPM) Indicator.
- Clear Acrylic with head puller and safety lid interlock to prevent cover opening during centrifugation.
- Inbuilt digital countdown timer 0-60 minutes .
- **Speed Capacity:** 6000RPM Max.
- **Speed accuracy:** +/-20
- **Time range:** 1~30 min.
- **Rotor capacity:** 6x15ml
- **Total volume:** 300ml
- **Head & Accessories :** 8x15ml angle head complete capacity with SS tube holder.
- **Power Supply:** 220-240 VAC, 50Hz.
- **Dimension :** 50 x 50 x 50 cm
- **Net. Weight :** 16 Kg (Approx.)

NOTE : DUE TO OUR CONTINUOUS R & D PROGRAM SPECIFICATIONS AND DESIGN MAY CHANGE WITHOUT PRIOR NOTICE.



C Cube Advanced Technologies
(An ISO-9001 & ISO-13485
Certified Company)

+91-80500 50512 +91-80-43021028

marketing@ccubeadvtech.com
ccubeadvtech@gmail.com

www.ccubeadvtech.com

Manufacturers, Suppliers & Dealers of Medical Equipments