

DENTAL CHAIR

Owing to its state of the art design, cushioned seat, user-friendly operations, and various tools that assist dentists in diagnosing and curing of a certain oral condition, this chair is widely in demand in the market.



Specifications

- Electrically operated dental chair with zero function.
- Chair mount unit 3 point modular system with hanging cords with ceramic spittoon.
- Chair mount unit with 1 airrottor, micro motor provision, three way syringe.
- Operating light with two intensity levels.
- Chair side spittoon and tubular water connection operated with feather touch switches.
- Common single foot control for both airrottor and micro motor.
- X-ray viewer and transparent water booster with regulator.
- Pneumatically operated high & low vacuum section.
- Joy stick foot control for chair up and down, back / front movements.
- Pneumatically operated doctor stool.
- All cushions like upholstery and operator stool made by PU moulded foam.

Operator Stool

- Pneumatic operating stool on five castor wheels
- Height adjustable between 500mm to 700mm
- Seat rotatable through 360°.

Compressor

- **Power :** ≥ 0.75 hp
- **Working pressure :** 8.00 Kg/c.m²
- **Displacement :** > 72 lit/Min
- **Operating Voltage :** 220 V + 10%
- **Tank Capacity :** 15 lit/50 lit
- **Noise Level :** <80dB

NOTE : DUE TO OUR CONTINUOUS R & D PROGRAM SPECIFICATIONS AND DESIGN MAY CHANGE WITHOUT PRIOR NOTICE.

