

Dermatology

Diode Laser



TABLE TOP MODEL



Vertical 808nm diode laser hair removal machine
TOWER MODEL

Specifications

Laser Type : Diode Laser

Laser wavelength : 808nm

Display : 10.4 inch color touch LCD screen

Spot Size : 12x12 mm²

Pulse width : 10-400ms adjustable

Energy : 10-200J / cm²

Frequency : 1-10Hz (continuous mode)

Cooling : water + air + semiconductor

Temperature of Probe : 0~4°C

Handpiece Cooling System : Sapphire TEC cooling+water cooling+Real-time monitoring

Water System Cooling : German high power condenser module+Fan cooling+Real-time monitoring

Water System Protection : Italian water pump, water flow, water level, water temperature etc.

Continuous Operation time : 24 hours

Laser bars : Germany imported 10 bars

Power : 600W

Size : 36x30x42cm

Voltage : AC220V±22V, 50/60Hz, 10A

Weight : 50KG

