

Dermatology

Elight

Crystal Erayzer - Multi Functional Laser Platform with IPL, Elight, RF, SHR, Q switch ND Yag Laser



Elight (IPL+ RF) System

Energy : 1-50J adjustable

Pulse : Single

Pulse Duration : 1-15ms

Pulse Repetition Rate : 1-10HZ

Handle : 1 Elight handle

Filters : 480/560/640nm standard (430nm,530nm, 590nm,690nm,750nm optional)

RF System

RF Energy : 0-50J adjustable

RF Frequency : 6MHZ

Handle : 1 RF handle with 3 probes

Spot Size : 20 / 28 / 35mm

ND YAG Laser

Laser Type : Q Switch ND YAG Laser

Wavelength : 1064nm,532nm

Infrared Aim Light : Infrared ray indicator

Energy : 0-1500MJ

Machine Specification

Power : 2000W

Cooling System : Wind+water+semiconductor+gel

Controller : 8 inch real color touch screen

Voltage : 220V(110V)/50HZ (60HZ)

Dimension of machine : 51x46x54cm

Dimension of package : 63x60x55cm

Applications

- Permanent hair removal
- Acne treatment
- Vascular Treatment
- Skin rejuvenation
- Wrinkle removal treatment
- Skin Tightening
- Tattoo removal
- Birth mark removal

