

# Dermatology

## Hand and Foot Unit

### Application

Psoriasis, Vitiligo, Atopic Dermatitis, Circumscribed Scleroderma and Other Skin Disorders.

### Salient features

- LCD JOULE Timer
- UVA & NBUBV treatment Selection provided in the controller.
- Large bright display are easy to monitor and for viewing from a distance.
- LCD display system provided to set the required joules
- Feather touch keys provided for easy operation
- LCD system provided to indicate the remaining time / joules to complete the therapy
- Provision for Usage hour monitoring in controller to enable timely replacement of tubes
- Four Panels with tilting upto 90°
- Unit Dimension : 12 x 85 x 165 cm excluding base
- Base dimension : 50(L) X 8(W) X 4(H)cms
- No. of packing boxes : 2Nos.
- Unit weight 80kg
- Choke : Philips Make
- Starters : Philips Make
- Tubes : Philips UVA/ Narrowband 2ft Tubes
- Aluminum reflector (100%) behind the tubes
- Main power supply :
- 230 V ± 10V
- Power cord : 3m, 6A
- Treatment covering area Dimensions 60x150 cm
- Safety Grill (MS) for UV tubes.
- Four 3 inch wheels to move the unit (2 with break 2w/o break)
- Handle provided for easily move the Unit to anywhere.
- Body of the entire unit :
- Made of Powder coated Mild steel
- Color : combination of 7035 (white) and DA gray (gray)
- Safety switch,
- Safety fuse and spare fuse
- Protection : safety earthing at 3places.
- Detachable power cord.
- UV Protective Goggle 2 nos, provided to avoid the UV rays.
- Tube to grill in between distance is 3 cm
- User manual provided to operate the unit easily



Variants	8 Tubes	12 tubes	16 tubes
UVA	HFU 8 UVA	HFU 12 UVA	HFU 16 UVA
NB UVB	HFU 8 NB	HFU 12 NB	HFU 16 NB
UVA + NBUBV	-----	-----	HFU 8+8