

# Dermatology

## Iontophoresis

It is a simple device very useful in palmer and planter as well as Ampits Hyperhidrosis. It provides for an effective procedure that yields excellent results relieving Hyperhidrosis of palms and soles and in ampits.

Plain Tap Water Iontophoresis in Palmo Plantar Hyperhidrosis

Iontophoresis is a treatment whereby ions of various soluble salts (e.g., Zinc, Chlorine, Iodine, Histamine) are introduced into the tissues by means of a constant electrical current - in brief, a form of electro-osmosis.

### Major Features

- An inexpensive machine, its operating costs are also very low.
- Its operation requires almost negligible prior training.
- Very effective, safe and gives superior results, thereby increasing the efficiency of the physician.
- Improves the quality of the medical care.

