


## Product Specification

<b>Product Name</b>	<b>Iontophoresis Unit</b>
<b>Model No.</b>	<b>CUT SWEAT</b>
<b>Make</b>	C Cube Advanced Technologies

Product	Product Description
	<p>Iontophoresis device is based on SMPS technology. It Works on Very low voltage &amp; current so it is a safest technology till date.</p> <p>It is very Affordable, safe to use, provide good results with economic cost.</p> <p>Ready to use and almost no installation procedure.</p> <p>Its operation requires almost negligible prior training. It is a Portable unit and fits to any smallest baggage very safe to use for hyperhidrosis.</p> <p>It is a simple device very useful in palmer and planter as well as Armpits Hyperhidrosis.</p> <p>Plain tap water can be used for Iontophoresis in palmo planter hyperhidrosis. Iontophoresis is a treatment where by ions of various soluble salts(e.g., Zinc, Chlorine, Iodine, Histamine) are introduced into the tissues by means of a constants electrical current - in breif, a form of electro-osmosis.</p>

Application	Iontophoresis is a treatment procedure that is used to treat hyperhidrosis disorder i.e. excess sweating in a person.	
Material	Unit is made up of ABS Plastic body	
Operation	Controller to set the treatment time and current flow in mA in Reverse/ Normal mode of operation.	
Product Dimensions	<b>Product Dimension</b> Product Dimension: 18 x8x6cm Product Weight: 0.950 Kg	<b>Dimension after Packing</b> Dimension: 30x35x8cm Weight : 1.2Kg
Power Source	230 V $\pm$ 10V, 50Hz	
Accessories	Pair of Patient Plate/ Earth Plate. Pair of Rollers or Underarm silicon flexible plates (Optional)at Extra cost Connecting cables Pair of Plastic Trays Bag to store and carry the unit User manual provided for easy operation of the unit	
Warranty	1 Year against Manufacturing defects	

### C Cube Advanced Technologies

Corporate Office: #406, 3rd floor, Malnad residency, 4th Main, Hanuma Hills Layout, Banashankri 3rd stage, Bangalore- 560 061  
 Factory: No.16/1, 9th Main Road, Opp to Taja Cake, Arehalli Village, Banashankari 3rd Stage, Subramanyapura Post, Bangalore-560061  
 Email: [info@ccubeadvtech.com](mailto:info@ccubeadvtech.com) / [marketing@ccubeadvtech.com](mailto:marketing@ccubeadvtech.com) / [ccubeadvtech@gmail.com](mailto:ccubeadvtech@gmail.com)