






Cosmetology

iTrans



-  40K Cavitation
-  Body RF Handle
-  Cryolipolysis (Medium) + Vacuum
-  Cryolipolysis (Big) + Vacuum
-  Cryolipolysis (Small) + Vacuum

The cryolipolysis handles can work at the same time. 3 handles for cryolipolysis

Technical Parameters

Display : 8 "TFT color touch screen

Voltage : 220V / 50Hz or 110V / 50Hz (according to the mark on the machine)

Power : ≤ 500w

Cooling liquid : Distilled water

Vacuum : 650mmHg

Pump flow rate : 60L / min

Vacuum pressure : 0-106Kpa

Freezing temperature : 0~, -10~, -15~

Fat exploding frequency : 35 - 40K

RF frequency : 5MHz

Packing size : 122cm × 59 cm × 51cm

Ambient temperature : 5 ~ + 40 ~

Relative Humidity : ≤ 80%



Accessories of this device

- Power supply line.....1pc
- Handle support.....1pc
- User manual.....1pc
- Cooling handle.....3pcs
- Cavitation handle.....1pc
- RF handle.....1pc
- Fuse replace.....2pcs
- Emergency switch.....1pc
- Anti freeze mask10pcs