

# Cosmetology

## Laser Comb



Product Name : Laser Comb

Model : REHAIR

Application : Hair Nourishment, Retention and Growth


Safe to use  
LLLT Technology

### A Revolutionary Treatment for Hair Loss & Hair Growth :

The Laser Comb provides low-level laser light energy that stimulates and enlivens hair follicles promoting new hair growth. With consistent use, most users start seeing benefits in 16 weeks. (Based on clinical studies)

#### Features :

- Portable, Handy with Unit Dimension: 21X5X4.5 cms
- 3 Laser Modules of 650nm with special Laser Multiplier
- Comb attachment for maximum laser penetration on the area to be treated
- 20 Minute Recommended Treatment Time
- 5A Adapter
- Very economically priced
- More manageability
- Increased Hair Growth
- Cessation of hair loss
- Faster growing hair
- Ease of use
- Operates at 230V  $\pm$  10V
- Designed as Box Piece for safe and dust free storage



## Rehair Plus



#### Features :

- Portable, Handy with Unit Dimension: 21X5X4.5 cms
- 10 Laser Modules of 650nm
- Comb attachment for maximum laser penetration on the area to be treated
- 15 Minute Recommended Treatment Time
- 5A Adapter
- Very economically priced
- More manageability
- Increased Hair Growth
- Cessation of hair loss
- Faster growing hair
- Ease of use
- Operates at 230V  $\pm$  10V
- Designed as Box Piece for safe and dust free storage