


PRODUCT SPECIFICATIONS

Product Name	HAIR GROWTH LASERS
Make	C Cube Advanced Technologies

Product	Model No	Product Description
	REHAIR	Unit Dimension: 21X5X4.5 cms After packing :24x11x9 cms Unit Package weight 1.0kg. Volumetric Weight: 1.4Kgs Adapter: 5 Volts 1Amps Laser Diodes: 3 Nos. 650nm 5mW Light output through special laser multiplier Comb Attachment: Fixed Main power supply: 230 V \pm 10V
	REHAIR PLUS	Unit Dimension: 21X5X4.5 cms After packing :24x11x9 cms Unit Package weight 1.0kg. Volumetric Weight: 1.4Kgs Adapter: 5 Volts 1Amps Laser Diodes: 10 Nos. 650nm 5mW Comb Attachment: Fixed Main power supply: 230 V \pm 10V
	LASER OPTIMUS	Unit Dimension: 29X23X20 cms After packing :34x27x25 cms Unit Package weight 4 Kg. Volumetric Weight: 6 Kgs Rechargeable Battery Operated with auto cut off feature after 21 minutes Laser Diodes: 28 Nos. 650nm 5mW Light output through special laser multiplier Cooling Fan provided for patient comfort Main power supply: 230 V \pm 10V Body of the entire unit: FRP

C Cube Advanced Technologies

Corporate Office: #406, 3rd floor, Malnad residency, 4th Main, Hanuma Hills Layout, Banashankri 3rd stage, Bangalore- 560 061
Factory: No.16/1, 9th Main Road, Opp to Taja Cake, Arehalli Village, Banashankari 3rd Stage, Subramanyapura Post, Bangalore-560061
Email: info@ccubeadvtech.com / marketing@ccubeadvtech.com / ccubeadvtech@gmail.com

	HAIR4SURE	Unit Dimension of laser hood: 38X34X19 cms After packing : 45X39X25 cms Unit Package weight of the hood 8 Kg and stand 15 Kg Volumetric Weight of the hood: 10 Kg and stand 30 Kg Timer provided with option to set the therapy time and once therapy is elapsed the unit will switch off automatically. Laser Diodes: 60 Nos. 650nm 5mW Light output through special laser multiplier Cooling Fan provided for patient comfort Main power supply: 230 V + 10V Body of the entire unit: Powder coated MS with FRP
---	------------------	--

Application	To promote natural permanent hair growth, Helps in longer retention of hair after PRP treatment, Mesotherapy, Hair Transplantation. Reduces Hair loss. Prevents hair thinning
Material	Powder coated Mild steel+FRP
Power Source	230 V \pm 10V, 50Hz
Warranty	1 Year against Manufacturing defects. Ware and tare parts are not covered under warranty. Battery warranty is for 6 months.

C Cube Advanced Technologies

Corporate Office: #406, 3rd floor, Malnad residency, 4th Main, Hanuma Hills Layout, Banashankri 3rd stage, Bangalore- 560 061
Factory: No.16/1, 9th Main Road, Opp to Taja Cake, Arehalli Village, Banashankari 3rd Stage, Subramanyapura Post, Bangalore-560061
Email: info@ccubeadvtech.com / marketing@ccubeadvtech.com / ccubeadvtech@gmail.com