

Cosmetology

Microdermabrasion



Beauty Queen tower model



Beauty Queen II

Model : Beauty Queen tower , Beauty Queen II

Category : Dermatology / Cosmetology

Application : Acne Removal, Scar Removal, Stretch Marks removal, Wrinkle removal & Skin Rejuvenation

Salient features

- Vacuum pressure upto 650mmHg
- Vacuum gauge provided to set the Vacuum Pressure and hour meter provided to know the usage hour.
- Unit Operation Mode: Diamond Tip and Crystal Powder (Beauty Queen Tower)
Diamond Tip (Beauty Queen II)
- Exhaust Fan provided to avoid the heat.
- Main power supply: 230V \pm 10V , 50 Hz
- Connecting Cable : 6A, 3CORE CABLE 3m long
- Unit Details:
 - ◆ The unit is made of MS. Drawer provided to keep the hand piece.
 - ◆ Color: combination of cream and blue. Fiber glass Finishing.
- User manual provided for easy operation of the unit.
- Warranty: One Year against Manufacturing Defects.

