

Cosmetology

MIFU - Slimming Machine



Specifications

Cryolipolysis + 40K Cavitation + Multi-polar RF + Lipo Laser

Technical Parameter

Voltage :	AC110V/220V 50-60Hz
Power :	< 400w
Temperature :	-14 °C~14°C
Pressure output :	0-900Kpa
Handle working size :	5.5*14cm
Screen :	8.4 inch touch screen
Body RF :	1Mhz
Lipo laser :	635nm-650nm
Cryolipolysis :	-14 °C~15°C
Cavitation frequency :	40Khz



Accessories List

Description	Qty.
Main machine (with cryolipolysis Handle)	1pc
Anti freeze mask	10pcs
Body RF	1pc
Cavitation handle	1pc
Laser pads	8pcs
Bandage	6pcs
Goggle	1pc
Water funnel	1pc
Power cord	1pc
Fuse	2pcs
Holder	2sets
Manual	1book