

NMI SYSTEMS



Needlefree Mesotherapy Injection System

Needle free Mesotherapy injection system directly acts on the skin to instantly enhance skin tissue penetration. During the process of electroporation, microspores of the skin gets activated and increases its absorption capacity, thereby enhancing the result of the induced medication.

Applications

- Freckle removal
- Skin whitening
- Skin toning
- Induce medicines, liquids with medicinal values or any beauty enhancement nutrition treatment
- Anti aging
- Skin rejuvenation
- Wrinkle treatment

Package

- Machine - 1pc
- Meso gun - 1pc
- Conductive probe - 1pc
- Serum bottle - 1pc
- The cover of serum bottle (White) - 1pc
- Power cord - 1pc
- User manual - 1pc


ACNE LIGHT ACNILL



The device adopts stainless steel head structure with heating feature along with blue light of 415nm known for antiseptic action and to eliminate the role of inflammation on the skin for the removal of facial acne. Device vibrates with 12000RPM to generate negative ions which play important role in repairing the skin.

Specifications

- **Power** : 2-AAA battery
- **Vibration** : 12000RPM
- **Material** : ABS + stainless steel
- **Net weight** : 34.5g
- **Blue wavelength** : 415 + 5nm , 415 - 5nm.



NOTE : DUE TO OUR CONTINUOUS R & D PROGRAM SPECIFICATIONS AND DESIGN MAY CHANGE WITHOUT PRIOR NOTICE.