







PRODUCT SPECIFICATIONS

Product Name	Phlebotomy Chair/ Blood Donor Chair
Product Description	Blood Sample collection Chair and Blood Donor Chair for Clinic and for Laboratories
Make	C Cube Advanced Technologies

Product	Model No.	Product Description
	P-200	Blood Sample Collection Chair Chair Dimension in inches: L: 26 D: 20 H: 18
	P-300	Blood Sample Collection Chair Chair Dimension in inches: L: 26 D: 20 H: 18
	P-500	Blood Sample Collection Chair with Material storage. Chair Dimension in inches: L: 36 D: 20 H: 18
	P-600	Blood Donor Chair with Manual back rest and foot rest Adjustment/Movement feature. This chair can be used from chair position to sleeping position Chair Dimension in inches: L: 70 D: 22 H: 26

C Cube Advanced Technologies

Corporate Office: #406, 3rd floor, Malnad residency, 4th Main, Hanuma Hills Layout, Banashankri 3rd stage, Bangalore- 560 061
 Factory: No.16/1, 9th Main Road, Opp to Taja Cake, Arehalli Village, Banashankari 3rd Stage, Subramanyapura Post, Bangalore-560061
 Email: info@ccubeadvtech.com / marketing@ccubeadvtech.com / ccubeadvtech@gmail.com

	P-700	<p>Blood Donor Chair with Manual back rest and Height adjustment feature, foot rest is fixed. Height adjustable rotatable armrests. This chair can be used from sitting position to sleeping position</p> <p>Chair Dimension in inches: L: 70 W: 34 H: 12 to 26 adjustable</p> <p>SPECIALY DESIGNED TO SUIT REQUIREMENTS OF BLOOD DONATION CAMPS (Can be easily folded and assembled in 5 to 10 min)</p>
	P-800	<p>Remote controlled electrically operated height adjustable sitting to sleeping position blood donor couch</p> <p>Height adjustable rotatable armrests</p> <p>Chair Dimension in inches: L: 70 W: 24 H: 22 to 34 adjustable</p> <p>Power source: 230V, 50Hz</p>
Application	Blood donor couch, Blood donor chair, Pathology labs, Blood donation camps	
Material	Powder coated Mild Steel	
Warranty	1 Year against Manufacturing defects. Ware and tare parts not covered in warranty.	

C Cube Advanced Technologies

Corporate Office: #406, 3rd floor, Malnad residency, 4th Main, Hanuma Hills Layout, Banashankri 3rd stage, Bangalore- 560 061
Factory: No.16/1, 9th Main Road, Opp to Taja Cake, Arehalli Village, Banashankari 3rd Stage, Subramanyapura Post, Bangalore-560061
Email: info@ccubeadvtech.com / marketing@ccubeadvtech.com / ccubeadvtech@gmail.com