

Dermatology

Portable UV Therapy unit

Application

Psoriasis, Vitiligo, Atopic Dermatitis, Circumscribed Scleroderma and other skin disorders.

Salient features

- Unit Dimension : 45 x 35 x 80 cm
- After packing : 48x40x90cm
- Unit weight : 10kg.
- Volumetric Weight after package : 12Kgs
- Choke : Philips
- tubes:
 - ◆ Philips
 - ◆ Tube size : 2 Feet.
 - ◆ Aluminum reflector (100%) behind the tubes
- Main power supply:
 - ◆ 230 V \pm 10V
 - ◆ Power cord : 3Yard, 6A

Varients	With Electromechanical timer	With LCD Digital timer
UVA Philips Tubes	PUVT UVA	PUVT UVA PLUS
NB UVB Philips Tubes	PUVT NB UVB	PUVT NB UVB PLUS



PUVT UVA/PUVT NB UVB with LCD Timer



PUVT UVA/PUVT NB UVB with Electromechanical Timer

- Optical Window Dimensions : 58 x 23cm
- Safety Grill for UV tubes.
- 2 inch nylon bushes provided at backside of the unit to use the unit in horizontal position.
- Handle provided for easy movement.
- Body of the entire unit :
 - ◆ Made of Powder coated MS.
 - ◆ Color : combination of 7035 (white) and DA gray (gray)
 - ◆ Switch provided for manual On / Off
 - ◆ Safety earthing at 2 points.
 - ◆ Detachable power cord.
 - ◆ Four 1 inch Bushes at the bottom of legs
- UV Protective Goggle provided to avoid the UV rays.
- Bag provided to avoid the dust.
- Tube to grill distance is 3cm.
- User manual provided to operate the unit easily