

PROCEDURE CHAIR HT-100 (3 actuator) / HT-200 (2 actuator)



Applications

General examination, dermatology and cosmetology procedures, Laser treatments, Hair Transplantation.

Specifications

- **Make & Model:** HT-100 (3 actuator) / HT-200 (2 actuator)
- **Material :** Epoxy Powder Coated Steel.
- Wired remote controlled Actuator based technology
- Body made of Powder coated mild steel
- Weight bearing capacity upto 150Kgs
- Height adjustable rotatable armrest
- Detachable handrests for sleeping position
- Tiltable height adjustable headrest
- PU foam for maximum patient comfort
- **Dimensions:** Length-6Ft, Width-2 Ft height adjustable from 24 inches to 38inches
- This model can be fitted with either 2 or 3 actuators to control the movements
- One of the actuator controls up and down movement
- Second actuator controls synchronous movement of back rest to foot rest.
- While the third actuators enable separate controls for backrest and footrest. The chair can be used for general examination, dermatological and cosmetologically procedures, Laser treatments, Hair Transplantation (**OPTIONAL**)
- **Product Dimension & Weight :** 72x24x24-38inch & 100Kg

NOTE : DUE TO OUR CONTINUOUS R & D PROGRAM SPECIFICATIONS AND DESIGN MAY CHANGE WITHOUT PRIOR NOTICE.



C Cube Advanced Technologies
(An ISO-9001 & ISO-13485
Certified Company)

+91-80500 50512 +91-80-43021028

marketing@ccubeadvtech.com
ccubeadvtech@gmail.com

www.ccubeadvtech.com

Manufacturers, Suppliers & Dealers of Medical Equipments