

Dermatology

Transkin – Skin Analyzer

- Nice appearance and innovative design. LED illuminator around lens.
- Imported lens with plated layer.
- 5.0 Mega pixels high resolution Korea sensor.
- Special DSP image processor, Optical Image Stabilizer.
- Adjustable focus to give clear image.
- Auto white balance and contrast adjustment, Color Temperature Filter.
- Dual image compare function.
- 3D-Negative capture model.
- Compatible with iris lens, hair lens and skin lens.
- Deliver clear and accurate images.
- Easy to operate.
- Maximum resolution: 2560x1920.
- OS : Windows XP, WIN2000, 2003, Vista, Win7, Win8.

Computerized
report generation
system

