

# Cosmetology

## Slimming Machine



### Technical Parameters

Power supply Input : 110V/220V 60HZ/50HZ

Power : 130W

#### 40K cavitation head :

Power supply output : 150V

Frequency : 40KHZ

Power : 25W

#### Vacuum+Bipolar RF head :

Power supply output: 66V

Frequency : 450K

Power : 95W



#### Body RF head :

Power supply output : 66V

Frequency : 550k

Power : 50W

#### Tripolar RF head :

Power supply output : 66V

Frequency : 450K

Power : 35W

#### Face RF head :

Power supply output : 66V

Frequency : 450K

Power : 35W

### Features & Functions

- The entire process is completed without the need of surgery and anesthesia.
- Does not cause the phenomenon of uneven skin.
- Does not cause the phenomenon of bleeding, swelling and blood stasis.
- Without side effects and the risk of gaining the weight back, the results are rather obvious.
- Non-invasive treatment does not affect the normal work and life.

Functions : Firms the arms, legs, thighs, buttocks, waist, back, abdominal muscles and reshape the body figure.

- Improve and treat the orange peel like fatty tissues of different serious levels With the effects of soothing, decompression and pain relief (such as joint pain, sciatica), promote the regular mechanical aerobic exercise of skin tissue, rapidly relieve stress and remove fatigue.
- Firm, lift and compact the skin and mould the face shape.
- Circulate and activate cells, promote metabolism and so that 100% more collagen and elastic fibers may be gained for improving the facial and body skin.

#### Package Includes:

- Multipolar RF head(Tripolar RF) X 1
- 40KHz Cavitation Head X 1
- Vacuum Bipolar Head x 1
- Bipolar RF Head for face X 2
- Power cord x 1

