

Cosmetology

Teeth Whitening Laser



Specifications

Dimension of frame : Height : 1420mm,
Arm length : 1250mm(max),
Wide : 630mm

Weight of machine : 18.5Kg

Input voltage : AC100~240V,50/60Hz

Output voltage : 4.7V

Power: 40VA

Wavelength of blue light : 460-490nm

Density of blue light : 200-300mw/cm²

Whiten time : 10-15minutes for one time,
two or three times for a treatment

size : 67*66*28CM

weight : 24KG



Advantages

- Our outstanding cool light technology is researched and developed by professional team, which is more effective than other teeth whitening system, ensure you get perfectly white teeth in short time.
- The specific wavelength accelerates the oxidization process and will quickly achieve whitening.
- The semiconductor light source will thoroughly prevent the disadvantage of halogen and harmful UV light and achieve truly luminescent whitening, the semiconductor has a guaranteed long life(over 50,000 hours).
- The LED head unit has a unique design, the light source has wide range and the light comes in at a perfect angle and radiates fully and uniformly.
- The distance between the LED head unit and cheek retractor will make sure the the LED head unit keeps clean, prevents cross infection, and meet high hygiene standards.
- All the parts are made of aluminum, which can last for longer use.

