

Cosmetology

Thermage



Specifications

Treatment Handle: 1 handle with 3 tips for body, face, eye respectively.

Type : Fractional RF

Heat : 40° - 60° centigrade

Voltage : AC110V / 220V, 50HZ / 60HZ

Power : 120VA

Fuse : 5A



Q Switch Nd Yag Tattoo Removal



Specifications

Laser type : Q switched ND:YAG laser .

Cooling : water cooling & air cooling

Spot size : 1-7mm(continuously adjustable)

Energy : <1000mj(multi pulse)

Wavelength : 1064nm, 532nm

Pulse width : 6-12ns

Input power : 500W

Gross weight : 30KG

