

Dermatology

UV Comb



Salient features of UV Comb

- UV Comb Unit Dimension : 5x 6 x30 cm.
- After packing : 10 x 38cm x 18cm.
- Unit Package weight : 1.0kg.
- Volumetric Weight : 1.4Kgs.
- Choke : 9 watts
- CFL Lamp :
 - ◆ Philips UV-A/NB UVB as per customer requirement
 - ◆ Power : 9 watt
 - ◆ 220- 240 V AC
 - ◆ Aluminum reflector (100%) backside of lamp.
- Main power supply : 230 V \pm 10V.
- Optical Window Dimensions : 10.5x 2.5 cm.
- Detachable Comb attachment.
- Body of the entire UV Therapy unit : Light Weight FRP
 - ◆ Color : combination of White and Gray.
 - ◆ On / Off Switch for Safety.
- UV Protective Goggle provided to avoid the UV rays.
- CFL UV Lamp to comb attachment tip distance is 3cm.
- 28 teeth in 2 rows of 14 each and teeth size is 3mm(Dia) x 17mm (Height).
- User manual provided to operate the unit easily. The unit can be used at home as per doctor's prescription. Hence the unit is also called as Mini Home UV Therapy Unit.

Salient features of UV Comb Plus

- Digital Timer with MM:SS (minutes:Seconds) format
- Large bright displays, easy to monitor and for viewing from a distance
- Digital display system provided to set the required time and to indicate the remaining time to automatically stop the unit
- Feather touch keys provided for easy operation.

UV Comb Plus



Application

Psoriasis, Vitiligo, Atopic Dermatitis, Circumscribed Scleroderma and Other Skin Disorders

