

## PORTABLE VETILATOR



**Standard Configuration :**  
Main Unit, Ventilation Hose, Patient Valve, Mask

Our Portable automatic oxygen respiration device (short-time ventilator) with additional inhalation facility. It can control and assist Patient whose weight is more than 15Kg, as well as non-invasive and invasive mechnaival ventilation.

### Application

Wild rescue, vehicle transport, long-distance transport, hospital transport.

### Features

- **Make & Model :** C Cube, AIRLOG3
- IPPV, Assist & Manual
- **TV range :** 75-2000ml
- **Breathing rate :** 5-40 BPM
- **Inspiration :** Expiration 1 : 1.67
- **Minute volume (MV) :** 3-20 LPM
- **Pressure limit :** 20-60 mbar, trigger at less than -2 mbar
- **PEEP range :** 5-20 mbars
- Alarm/warning once patient stop breathing
- Gross weight not over 1.3kg (just only the machine, not include tube and etc.)
- **IP class :** IPX4
- **Battery working time :** 10 hours



## MEDICAL VETILATOR



**Standard Configuration :**  
Ventilator, Oxygen Pressure Reducer, High-Pressure oxygen hose, Lung Simulator, Breath Circuit, Mask, Power Adapter.

Our Medical Ventilator is used for First-Aid, Ambulance, Emergency situation and patients' transport in hospital.

It can be used for child and adult patients.

### Features

- **Make & Model :** C Cube, MV-100
- Multiple functions, easy to take.
- The main parts adopt good-quality components, accuracy and reliable.
- Safety protection function.
- **Environmental Requirements:**
  - + **Relative humidity:** 15%~95%
  - + **Ambient temperature:** -18~ 50°C
  - + **Atmospheric pressure:** 70~106kPa
  - + **Gas source:** 0.28~0.6Mpa 60l/min
  - + **Power:** AC100V-240V
  - + **Li-battery:** 12 V, rechargeable
- **Classification:**
  - + Class, Type B—by protection against electric shock.
  - + IPX4—by protection against harmful ingress of water, particulate matter.

NOTE : DUE TO OUR CONTINUOUS R & D PROGRAM SPECIFICATIONS AND DESIGN MAY CHANGE WITHOUT PRIOR NOTICE.



**C Cube Advanced Technologies**  
(An ISO-9001 & ISO-13485  
Certified Company)



+91-80500 50512



+91-80-43021028



marketing@ccubeadvtech.com  
ccubeadvtech@gmail.com



www.ccubeadvtech.com

Manufacturers, Suppliers & Dealers of Medical Equipments